## GALEAS™ Hereditary Plus

A clinically validated NGS panel with optimized bioinformatics for analyzing germline mutations associated with hereditary cancers

#### **Highlights**

#### Enhanced clinically relevant content for assessment of hereditary cancer variants

Expertly curated coverage for the key clinically relevant regions of 146 genes associated with a predisposition for hereditary cancer including breast, prostate, Lynch syndrome and Wilms tumor.

#### A single workflow for all clinically relevant variant types, sample types and cancers.

Validate and run one workflow for all hereditary cancers, regardless of sample type (e.g. blood or saliva) and confidently call all variants including a wide range of CNVs in genes such as APC, MHS2, BRCA1 and PMS2 without the need for additional MLPA.

#### Supported by GALEAS Analysis software.

Optimised for the GALEAS panels, our cloud-based bioinformatics pipelines deliver accurate calling across all variant types associated with inherited cancers.

#### Introduction

Between 5-10% of all cancers, including cancers of the breast, ovary, uterus, prostate, and gastrointestinal system can be accounted for by hereditary cancers. The identification of individuals who are at an increased risk of developing inherited cancer is dependent upon the ability to accurately identify the genetic variants associated with heritable cancer syndromes. The ability to perform a comprehensive evaluation of the germline variants is key to understanding their association with cancer predisposition. This can provide a cancer risk assessment, and guide the implementation of additional screening and surveillance which may in turn result in an early diagnosis and guide treatment for both themselves and their families.

### Next Generation Sequencing (NGS)-based multigene panel for comprehensive profiling of heritable cancers.

The use of targeted NGS multigene panels to provide a comprehensive analysis of cancer susceptible genes has proven to be a clinically viable option for many laboratories. It allows researchers to profile known genetic associations for hereditary cancer regardless of sample type (blood or saliva) or cancer type.

However, many NGS panels struggle to identify key hereditary cancer copy number variants (CNVs), such as single exon BRCA1 or BRCA2 alterations or those CNVs involved in Lynch syndrome, and require additional multiplex ligation dependant probe amplification (MLPA) analysis to detect them.

# GALEAS Hereditary Plus panel design

GALEAS Hereditary Plus has been designed to target germline mutations in 146 genes associated with an increased risk of developing hereditary cancer. These genes have been selected to cover not only the common hereditary cancers like breast or prostate, but also the rarer hereditary cancer types like Phaeochromocytoma and paediatric cancers such as Wilms tumor.

Table 1: Genes included in key guidelines associated with risk of developing hereditary cancers and included in the GALEAS Hereditary Plus panel (see appendix for full gene list \* Wilms tumor only).

, , , , , , , , , , , , , , , , , , , ,			
Cancer type	Recommended genes for screening included in GALEAS Hereditary <sup>Plus</sup>		
Breast	ATM, BARDI, BRCA1, BRCA2, CDH1, CHEK2, NBN, NF1, PALB2, PTEN, STK11, TP53		
Colon	TAPC, AXIN2, BMPRIA, CHEK2, EPCAM, GREMI, MLHI, MSH2, MSH6, PMS2, MSH3, MUTYH, NTLHI, POLDI, POLE, PTEN, RNF43, SMAD4, STKII, TP53		
Renal	BAP1, FH, FLCN, MET, SDHB, VHL		
Ovarian	ATM, BARDI, BRCA1, BRCA2, CDH1, CHEK2, NBN, NF1, PALB2, PTEN, SKT11, TP53, RAD51C, RAD51D		
Prostate	ATM, BRCA1, BRCA2, CHEK2, MLH1, MSH2, MSH6, PALB2		
Gastric/GIST	CDH1, KIT, PDGFRA, SDHC, SDHD, SDHA		
Brain	APC, ATM, MLH1, MSH2, MSH6, PMS2, TP53		
Sarcoma	EXT1, EXT2, MTAP, NF1, RECQL4, SQSTM1, TP53		
Paediatric*	CDKNIC, CTR9, REST, TRIM28, WTI		

The design has been carefully curated to ensure that all clinically relevant exons are covered, including selected non-coding regions such as BRCA1/2 5' UTRs and the APC promoter. The GALEAS Hereditary Plus also includes a panel of tracking SNPs for patient identification purposes. When combined with the GALEAS analysis software this panel allows the sensitive and specific detection of SNVs, INDELs and CNVs.

Table 2: GALEAS Hereditary Plus panel specifications (\* for gDNA only)

Parameters	Specification	
Enrichment method	Hybridization and Capture	
Number of genes	146	
Capture Panel size	809 Kb	
Sequencing platform	Illumina	
Targets	Genes associated with hereditary cancer	
Variant types	SNVs, CNVs and INDELs	
Input DNA requirements*	10-200ng	
Sample type	gDNA from blood or saliva	
Multiplexing guidance for sequencing*	1 Million reads per sample required to achieve 100x. This equates to 0.2Gb per sample.	

## GALEAS Hereditary Plus panel performance

### Superior precision and recall ensure confident calling of SNV and indel variants.

GALEAS Hereditary Plus was validated across 437 SNPs using commercially available reference control NA24385.

Table 3: GALEAS Hereditary Plus SNV and indel recall across 4 replicates of reference standard NA24385.

	Recall
SNV	99.78%
Indel	100%

### Superior precision and recall ensure confident calling of copy number variants

To evaluate the sensitivity of CNV genotyping with GALEAS Hereditary Plus, the panel was run using NIBSC Lynch Syndrome MLPA cell lines. All CNVs were detected with 100% recall and precision when using sex matched control pools (Table 4).

Table 4: Recall and precision statistics for copy number alterations (CNVs) in NIBSC reference controls using the GALEAS Hereditary Plus panel.

CNV	Genotypic sex	CNV type	Detected
Copy normal	Male	Copy neutral	Yes
MSH2 deletion exons 1-6, heterozygous	Male	Multi-exon deletion	Yes
MSH2 deletion exon 7, heterozygous	1 Male	Single exon deletion	Yes
MSH2 deletion exons 1-2, heterozygous	female	Multi-exon deletion	Yes
MSH2 deletion, exon 1, heterozygous	Male	Single exon deletion	Yes
MLH1 exon 13 amplification (3 or more copies)	Female	Multi-exon amplification	Yes

## Panel performance specifications

GALEAS Hereditary Plus panel design delivers a high percentage of on-target reads, lower duplication rates and more consistent vertical coverage with 99% of targets covered at 30x or more (Table 5). This exceptional technical performance delivers high recall and precision across more variants associated with hereditary cancer, including CNVs, than a leading competitor without significantly increasing sequencing costs.

Table 5: Sequencing metrics for GALEAS Hereditary Plus. Compared with another commercial alternative, GALEAS Hereditary Plus delivers 100% more content (including CNV probes) for less than 10% more sequencing.

Key Quality Indicator	GALEAS Hereditary <sup>Plus</sup>	Company I
Number of genes	146	113
Capture panel size (kb)	809 kb	403 kb
GB required for mean 100x coverage	0.2 Gb	0.12 Gb
Percentage coverage >30x	99%	96%
Percentage on or near bait	81%	61.51%
Percent duplication	2.0%	8.99%
SNV recall	99.7%	98.1%
INDEL recall	100%	97.2%

## GALEAS Hereditary Plus clinical validation

The clinical utility of GALEAS Hereditary Plus was assessed using research collaborator samples, including gDNA from 64 patient blood samples with orthogonal data.

#### SNV recall and precision

SNV recall on clinical samples was shown to be 100%, across a wide range of alteration types, including small and large (>10bp) INDELs.


Table 6: SNV recall on clinical samples for GALEAS Hereditary Plus

ID	Gene	HGVS coding	HGVS protein	Genomic position
22	BRCA1	c.1175_1214del	p.Leu392fs*5	chrl7:43094317
23	BRCA1	c.1175_1214del	p.Leu392fs*5	chrl7:43094317
64	MSH2	c.942+3A>T	P.?	chr2:47414421
65	PMS2	c.736_741delins TGTGTGTGAAG	p.(Pro246Cysfs*3)	chr7:5997389
66	MLH1	c.1946dupC	p.(Leu650Phefs*14)	chr3:37048561
67	MSH2	c.1213_1217dup	p.(Leu407Thrfs*7)	chr2:47429877
68	MSH6	c.3562_3563del	l p.(Serll88Tyrfs*5)	chr2:47805623

#### Copy Number Variants (CNVs) and pseudogenes

GALEAS Hereditary Plus accurately identifies CNVs from single exons to whole genes in key cancer syndrome susceptibility genes. When combined with the GALEAS Analysis Software, the GALEAS Hereditary Plus panel provides:

- An analytical sensitivity of 100%
- An analytical specificity of 93.5%
- The ability to detect mosaic copy number variations in key genes such as APC and TSC2 (see Figure 1, E and F)
- Capability to accurately distinguish between PMS2 and PMS2CL pseudogene (see Figure 1, C)



### **GALEAS** analysis software

GALEAS Analysis Software is a cloud-based set of optimised bioinformatics pipelines which provides accurate calling of SNVs, INDELs and a wide range of CNVs from single exons to whole genes. The pipeline can distinguish the PMS2 and PMS2CL pseudogenes and provides the ability to detect mosaic copy number variation in key genes such as APC and TSC2.

In addition, the GALEAS Analysis Software provides an easy-to-use method of uploading batches of FASTQ files and downloading the results, with just a few clicks.


#### **Panel of Normals**

To minimize background noise, improve CNV calling and reduce costs, the GALEAS Analysis Software leverages an in-built 'panel of normals'. Using a large cohort of clinically normal samples, processed with the GALEAS Hereditary Plus library prep and sequenced using an Illumina sequencing workflow, this data set provides a baseline to call CNVs, dramatically improving the accuracy of CNV calling.

# Streamlined workflows; quick and easy protocols

GALEAS Hereditary Plus enables laboratories to validate and run a single comprehensive workflow to profile all hereditary cancer types, reducing turnaround time, validation and operating costs. This can significantly reduce laboratory validation and operating costs. Validate and run one workflow for all hereditary cancers and confidentially call all variants including a wide range of CNVs in genes like APC, MHS2, BRCA1 and PMS2; potentially eliminating the need for MLPA for clinical identification and reporting of CNVs.

The GALEAS Hereditary Plus workflow is simple and easy. Taking less than 10 hours, with less than 2 hours handson time, it is designed with multiple stop points to provide flexibility within laboratory processing. Library preparation can be run manually or automated up to 96 samples in a single run. Indexes are available for up to 384 samples to allow for flexible batch sizes and scalability across all Illumina benchtop sequencers.



### **Summary**

GALEAS Hereditary Plus provides an expertly curated, comprehensive NGS solution for the analysis of genes previously linked with cancer predisposition syndromes. The enhanced probe design, comprehensive coverage, high coverage and uniformity allows the accurate and sensitive detection of SNVs/INDELs and CNVs. Combining this with the GALEAS Analysis software provides a simple and easy sample to analysis workflow. GALEAS Hereditary Plus provides a highly efficient, targeted sequencing and analysis solution to allow the detection of variants associated with cancer predisposition syndromes.

#### Learn more

To learn more about GALEAS Hereditary Plus and to download the protocols, application notes and white papers please visit: **www.nonacus.com**.

#### References

1. Ngeow, J., Eng, C. Precision medicine inheritable cancer: when somatic tumour testing and germline mutations meet. npj Genomic Med 1, 15006 (2016).

#### Ordering information

#### **Product**

GALEAS Hereditary Plus 16 samples GALEAS Hereditary Plus 96 samples

#### Catalogue No.

NGS\_GAL\_HCP\_FR\_16 NGS\_GAL\_HCP\_FR\_96

#### **Nonacus Limited**

Quinton Business Park 11 Ridgeway Birmingham B32 IAF

info@nonacus.com

© 2023 Nonacus Limited. For research use only. Not for use in diagnostic procedures